Introduction

Fluorescence

- Three-step process:
- 1) Excitation by absorbing radiation (λ_{EX})
- 2) Vibrational relaxation and internal conversion to lowest vibrational level of the excited state
- 3) Emission of a photon of lower energy (λ_{FL})

RAMAPO

COLLEGE

Fluorescence Spectroscopy

- Polychromatic light is wavelength-selected by a monochromator and used to excite the fluorophore
- Fluorescence emission is collected at 90° relative to the excitation and directed towards an emission monochromator and a detector

Coumarin-120 (C-120) Fluorescent Probe


- Hydrophilic molecule
- Intensity depends on environmental conditions including pH, solvent polarity, and viscosity
- Useful in cell tracking applications


Reverse Micelles


- Surfactants are composed of a hydrophilic head and hydrophobic tail(s)
- Aggregate into ordered structures, such as micelles, due to entropic processes
- Critical Micelle Concentration (CMC) lowest concentration of surfactant required for aggregation
- Reverse micelles are composed of hydrophilic heads with hydrophobic tails extending outwards
- Applications include targeted drug delivery, and protein isolation and refolding
- Hydrophilic fluorophores may be encapsulated in the core and used to determine the CMC by measuring fluorescence intensity


Aerosol-OT (AOT)


- Most widely studied surfactant molecule
- Capable of forming microemulsions
- Excellent in encapsulating hydrophilic molecules





















Determination of Critical Micelle Concentration of Aerosol-OT in a Two-Solvent System by Fluorescence Spectroscopy

Alexa Hacopian and Dr. Loraine Tan

School of Theoretical and Applied Science, Ramapo College of New Jersey, Mahwah, NJ


Experimental

Spectrofluorometer Parameters

Results and Discussion

Determination of CMC

- CMC study



Conclusions & Future Work

- resolution of spectra
- deviation

Acknowledgements

- Ramapo College of NJ



 λ increments: 1 nm • Excitation bandpass: 3 nm Emission bandpass: 3 nm

• Literature CMC of 1.09 mM was used as a center point for

• Experimental CMC was determined to be 1.14 mM (±0.35mM)

CMC of AOT in Hexane

• Concentration of C-120 should be reoptimized to improve

• Deposition method using THF was successful in encapsulating the fluorophore, should be optimized to lower standard

Another fluorophore, soluble in a nonpolar environment, such as dansyl chloride, can be used to determine CMC

References

 Olesik, Susan; Miller, Curtis J. Critical Micelle Concentration of AOT in Supercritical Alkanes; *Langmuir* 1990, **6**, 183-187 • Wagner, Brian. The Use of Coumarins as Environmentally-Sensitive Fluorescent Probes of Heterogeneous Inclusion Systems; *Molecules*, 1990, **14**, 210-237. • Telmo, J.V. Prazeres, Mariana Beija, Fabio V. Fernandes, Paulo G.A. Marcelino, Jose Paulo S. Farinha, J.M.G. Martinho. "Determination of Critical Micelle Concentration of Surfactants and Amphiphilic Block Copolymers Using Coumarin-153". Inorganica Chimica Acta 381 (2012) 181-187

Ramapo College of NJ, TAS Research Honors Program Aleksandar Goranov, Dalton Kaye, Michael Knef, and Brian Gural, Tan Research Lab Alumni (2015 - 2017) • Dr. Robert Mentore, Carol Ichinco, and Elizabeth Miloscia,